



A new age in packaging products presents a whole new area of sealing problems; specialised packaging materials for the pharmaceutical industry; all manner of blister packs for all manner of products; improved food and hygiene standards for flexible packaging materials make testing more necessary for the small confectionery company right up to the major pharmaceutical producer.

Options available:

- Heated upper jaw, silicone rubber lower jaw (supplied as standard).
- Heated upper and lower jaws (flat faced or crimped).
- Teflon coating.
- Blister pack sealing platen and tray holder.
- Gradient sealing with 3 or 5 jaws.
- Split sealing head and control.
- Alarm package for temperature and pressure.

Specification:

- Temperature:** Via self tuning PID control 0-400°C, 1°C resolution, type J thermocouple.
Accuracy +/- 0.5%.
- Pressure control:** Via precision regulator controller, 0-100 PSI (other units) repeatability 0.02 PSI indication via digital meter and strain gauge transducer.
Accuracy +/- 0.25%.
Operating accuracy +/- 5%.
- Dwell time:** Via digital quartz timer, range 0.1 to 999.9 seconds.
Accuracy +/- 0.01%
Operating accuracy +/- 5%.
- Jaws:** 300 x 25 as standard but can fit 350 x 25mm.
- Supply:** 220/240 VAC single phase
50/60 Hz (110/120 VAC available on request).



BS EN ISO 9001 Certificate No. 6220

Unit 46 Golds Nursery Business Park, Jenkins Drive, Elsenham,
 Bishop's Stortford, Herts CM22 6JX.
 Telephone: +44 (0) 1279 817171 Fax: +44 (0) 1279 815743
 e-mail: rdmtest@aol.com Internet: www.worldsites.net/rdm